


SAMBA HAIR

Hair Shine Analysis for Claims & Research



Polarization Imaging
Shine / Chroma / Diffused




Measurement
Automation

Claims Compliant
Hair Shine metrics



Increase in hair shine may be one of the most sought-after properties in the hair care industry, whether it be for shampoo, conditioner, oils or hair sprays. However, correctly quantifying such a subjective parameter is not easy.



The **SAMBA HAIR System** uses polarization-imaging to separate the different components of light that are reflected from hair swatches and extract **valuable data for research and claims of luster properties of hair**. The system includes an imaging sensor head and an in-vitro enclosure. The hair samples are positioned on a cylindrical mount that automatically moves in and out of the instrument.

SPECIFICATIONS

VISION	12-bits depth 3840 x 2160 px ² (camera) 1570 x 280 px ² (single swatch) White LED Illumination bar
HAIR SWATCH	Any color 0.5" (12.5 mm) wide 7" (180 mm) long
SOFTWARE	SAMBA HAIR 4.0 or higher
ACQUISITION	6 swatches max per cylinder 10 positions / measurement per swatch 5 to 20 minutes per full cylinder Shine / Chroma / Diffuse images
DATA	Luster coefficients (BNT, Reich-Robbins) Shine / Chroma / Diffuse Images and Angular distribution Profiles parameters (Intensity, band width, overlapping, etc.) Excel Export
SYSTEM	Size : 24" x 19" x 10" (61 cm x 49 cm x 26 cm) 110/240 VAC 50/60 Hz


SAMBA HAIR SOFTWARE

Using our polarization imaging system, the SAMBA Hair software can reconstruct the specular, shine, chroma and diffusion images of your hair swatches.


After processing, luster coefficients used for claims in the cosmetics industry as well as more extensive data (angular distribution of intensity, shine band parameters, etc.) are provided.



Specular and Diffusion images separation



Angular distribution
Shine / Chroma / Diffusion



With the SAMBA HAIR Software, you can easily :

- Acquire your images with a simple organization system for your studies
- Get luster coefficients widely used for claims in the cosmetic industry
- Compare images and numerical data of multiple samples
- Export all your data including images, graphs and tables